

# Comparing Characterization Methods for Investigation High-Temperature Bending Behavior of Melted Thermoplastic Prepreg using DMA and Rheometer

Luohaoran Wang<sup>1</sup>; Griffin Chase<sup>1</sup>; Logan Oster<sup>1</sup>; Jon Abramson<sup>2</sup>; Malcolm Lane<sup>3</sup>; Alireza Forghani<sup>3</sup>; Marco Salviato<sup>4</sup>; Navid Zobeiry<sup>1\*</sup>

1 University of Washington, Materials Science & Engineering, Seattle, WA, U.S.A.

2 Convergent Manufacturing Technologies US Inc., Seattle, WA, U.S.A.

3 Convergent Manufacturing Technologies Inc., Vancouver, BC, Canada.

4 Aeronautics & Astronautics Department, University of Washington, Seattle, WA, U.S.A.

American Society for Composites (ASC)  
38<sup>th</sup> Annual Technical Conference



***BE BOUNDLESS***







# Background

- Prepreg bending is key deformation mechanism during automated fiber placement (AFP) and hot drape forming
- Prepreg's bending response strongly depend on
  - Temperature
  - Loading rate





IMAGE CREDIT TO DLR GERMAN AEROSPACE CENTER



Farnand et al. 2017

# Convergent's Pure Bending Fixture

Relative work have done with rheometer by Convergent Manufacturing Technologies for thermoset prepreg

## Key features:


- Double roller design:
  - Minimizing frictional and membrane effects
  - Providing pure bending moment
- Fixed sample length
- Vertical configuration
  - To eliminate the effect of gravity
- Large to small bending radii
  - To reproduce realistic processing conditions

## Fully definable test conditions:


- Temperature
- Bending Rate
- Bending Angle and Radius of curvature
- Loading/Relaxation/Unloading

### Pure Bending fixture

inspired by Sachs (2014)





### Top View schematic of fixture



*\*Lane et al. 2023*

# Convergent test data for a thermoset prepreg




- Loading stiffness decrease with temperature increase
- Relaxation happens earlier with temperature increase





\*Lane et al. 2023

# Convergent's Bending Constitutive Model

- Fiber micro buckling is observed in bent samples
- This inspired a non-linear viscoelastic model that is represented as a hinge on a viscoelastic bed



Fibres buckle on the inner surface to accommodate the deformation.






# DMA Bending Method Development

## Dynamic Mechanical Analysis (DMA)

- built-in bending test capability and control
  - Temperature
  - Strain Rate
- Measure viscoelastic properties
  - Dynamic test
  - Relaxation test
  - Creep test
- In this research, we use DMA TA 850 to develop a new test method




# Material




Property of Toray TC1225 LMPAEK	Standard Value
Fiber areal weight (FAW)	145 g/m <sup>2</sup> (4.28 oz/yd <sup>2</sup> )
Resin content by weight (RC)	34%
Consolidated ply thickness (CPT)	0.14mm (0.0054 in.)
Glass Transition Temperature (T <sub>g</sub> )	147 °C (296.6 °F)
Melting Temperature (T <sub>m</sub> )	305 °C (581 °F)

## Transition Temperature



## Nominal bending behavior







# Test Methods

## 3-point bending/ 4-point bending schematic with moment & shear force distribution

- Non-pure bending moment with shear force
- Plastic hinge caused at contact point





- Isolate pure bending moment
- Constant curvature





# DMA Tests

- 3-point bending outer fixture
  - roller outer fixture could not support
  - customized clamp outer fixture with no span length change
- Developed customized 4-point bending fixture
  - To better isolate pure bending moment



# Final Fixtures

- Custom 4-point bending fixture
- Dual flat clamp outer fixture with Kapton tape
- Release agent applied on contact points



# Bending Calculation

According to bending theory, while neglected gravity, the beam deflection can be described by

$$EI \frac{d^2y}{dx^2} = M = \frac{F}{2} * x - \frac{F}{2} * \left(x - \frac{L}{4}\right) - \frac{F}{2} * \left(x - \frac{3L}{4}\right) \quad (1)$$


As a viscoelastic material, Pseudo viscoelastic *cure hardening instantaneously linear elastic* (CHILE) approach can be applied and the initial Euler-Bernoulli beam equation can be transformed to

$$dM = E'(\theta)I \times d\kappa \quad (2)$$

Final equations at  $\frac{L}{4}$  can be reached to approximate the melt bending behavior of the material :

$$M = \frac{FL}{8} \quad (3)$$

$$\kappa = \frac{12y}{L^2} \quad (4)$$




Where R is reaction force;  
F is load; G is gravity;  
E is Modulus;  
L is span length of specimen;  
b is width of specimen;  
E' is instantaneous elastic modulus;  
I is second moment of area;  
y is deflection;  
x is assumed distance




# DMA Results

- Data has been processed via the previously shown method
- Minimal change in bending rigidity up to about 300°C
- Notable stiffness drop at 305°C
- A softening response at 305°C
- Significant stiffness drop at 310°C





# Data Comparison

- Pure bending fixture in Rheometer allows tighter radius for testing
- Both methods indicate melting point around 305°C to 310°C
- Similar temperature-dependent rigidity behavior (Relaxation not include) is shown by both methods
- Attributed to buckling of fibers on compressive side of sample



# Material Behavior and Limitation of DMA

- During the bending process, micro-buckling can happen before relaxation (aligns with previous Convergent 's work)
  - Higher displacement
  - Higher applied load
- Curvature is limited to around  $20 \text{ m}^{-1}$



# Conclusions



- The developed 4-point bending fixture successfully isolates pure bending better than 3-point bending fixture
- Despite differences in loading application, sample geometry, and data reduction approach, similar results were obtained using both the Rheometer and DMA
- DMA test can characterize rigidity behavior
- DMA test has a minimum test radius around 5cm




# Acknowledgement

---

We wish to express our deepest gratitude to our sponsors, NASA, the Air Force Research Laboratory (AFRL), and Convergent Manufacturing Technologies US for their financial support.

Our sincere thanks also go to Convergent Manufacturing Technologies Canada for their technical assistance and collaboration.

Additionally, we extend our appreciation to Toray Advanced Composites for their generous material donation.



UNIVERSITY of WASHINGTON

**THANKS FOR COMING!**  
**Any Questions?**

